

# **Hohenstein**

# **Quality Label "Fit"**

## Objective

With the help of the independent product endorsement "Best suited to the material and functional fit", you can guarantee a consistent fit for your customers – right across all your collections.

## Description

- The Hohenstein quality label is particularly suitable for NOS goods.
- The first size set to be examined will be transferred from the manufacturer to Hohenstein to prove fit.
- The fit is tested both in new condition and after care from selected sizes and one sample washed as described on the care label.
- The assessment is made by the department fitting, development and pattern, using fitting models whose body measurements corresponds to "Size Germany" survey or the body-measurements table from the customer .
- If the size set requires the "pass", the product can be awarded with the Hohenstein quality label at the POS.



## Requirements for test samples

A product that reaches the highest assessment level both before and after washing can be awarded the Hohenstein quality label “fit”. This assures that the product fulfills the following requirements:

- The garment must fit with regard to all length and width measurements
- Freedom of movement must be guaranteed in accordance with its purpose and function
- The fit and appearance must not be impaired by defects in workmanship
- Children’s clothing must conform with DIN EN 14 682

## The scope of service

- Label printing template valid for 1 year
- Test report incl. assessment
- On request, a clarifying meeting including fitting/telefitting

## Customer benefit

- The neutral assurance of product quality by an independent institution can be shown to the end user at the POS in the form of a hangtag
- Quality assessment in new condition and after care
- Effective prevention of complaints at the POS

## What we need from you

- One size set in the original material incl. care instructions

## Label extension

- To extend the label for another year, our inspector picks up the latest size set for re-testing before the expiry of the period of validity

## Time required

- regular service: 10 working days
- express service: 5 working days